



Powdershape & Fibreshape



연구용 및 실시간 공정형 온-라인 분체 파우더 및 파이버 입자 입도 형상 영상 분석기



특징 및 장점

1. 스위스 IST AG사 섬유소 형상 영상장치 Fibreshape systems

Fibreshape 은 여러가지 섬유소 파이버 재료의 형상 분석 및 품질관리 (Quality Control: QC)를 위해 편리하게 사용할 수 있는 장치입니다. 파이버셰이프 기능, 사용자 인터페이스 및 신속한 분석 기능에 있어서 세계 유일의 Unique 한 모델이며 샘플 시료 조제에서부터 재료 인증서 발행까지 다양한 용도로 사용될 수 있습니다.

The measurement range is from **2 um to 30 cm fibre length**.

Fibreshape - 섬유소 시료를 손으로 스캐너에 조금씩 투입하는 수동 모델

Fibreshape automatic - 특허받은 특수 피더 장치를 이용하여 섬유소 시료를 자동으로 스캐너에 조금씩 투입시키는 전자동 모델 .

2. Two kind of Powdershape systems exist

Powdershape - where the samples are distributed stepwise by hand on the scanner

Powdershape automatic - where the samples are feded automatically by the sample feeder on the scanner. (sample feeder: patent pending).

Powdershape is a characterisation system for powders and generally all kinds of particles especially developed for quality inspection. It is unique in its set of features, simple user interface, and rapid characterisation process from the sample preparation to the printed certificate. The measurement range is from **1 um to 20 cm particle size**.

3. A new simply product "Microshape" has been issued for powder characterization that has just two measuring ranges

2 - 150 um using the digital microscope with 14600 dpi

10 - 1000 um using an 3600 dpi film scanner

연구용 및 실시간 공정형 온-라인 분체 파우더 및 파이버 입자 입도 형상 영상 분석기

Fibreshape

1. The Fibreshape system

Non-wovens (industrial fibres) control and characterisation system for fibre applications. It is unique in its set of features, simple user interface, and rapid characterisation process, from the sample preparation to the printed certificate.

The measurement range is from 2 μm to 30 cm fibre length.

Two kind of Fibreshape systems exist:

- **Fibreshape** - where the samples are distributed stepwise by hand on the scanner
- **Fibreshape automatic** - where the samples are feeded automatically by the sample feeder on the scanner. (sample feeder: patent pending).



The very simple handling of "Fibreshape automatic"

1. Simply distribute your fibre samples on the board
2. Enter a description and select fibre length, fibre thickness or other reports with a single click!
3. All the specimen are scanned and characerized automatically. Obtain the final reports or insert results into quality chart tabase

Your advantages:

- A very quick measurement method
- Easy and fast specimen/sample preparation
- Minimisation of possible operator errors
- Consistent and reproducible results
- Very attractive price/performance relation
- Reliable and rapid industrial method for quality inspection


2. Fibreshape features

Key features

- Quick measurement of fibre size and shape distribution
- Length measurement of short fibres (flock, linters and other fibres with a length/thickness ratio < 200)
- Fibre orientation measurement
- Characterisation of optical properties: colour and transparency
- Measurement of trash/dust content (consider mixtures for more details)
- Processing and storage of test data (for further external use and analysis)
- Integration for SAP or similar systems available
- Reliable and rapid industrial method for quality inspection

Fibres of a fineness range from 5 μm up to 5000 μm (5 mm) can be measured.

Simple user interface



To facilitate the fibre characterisation, the user interface exposes only 3 steps:

1. Image acquisition e recognition and measurement
3. Analysis by means of statistical evaluation

Technical data

The Fibreshape System consists of:

- The Fibreshape software
- State-of-the-art image analysis computer running Microsoft Windows NT/2000/XP
- A slide and/or flatbed scanner certified for Fibreshape
- Standard office suite

3. Fibreshape applications

- Wool
- Cotton
- Chipboard and paper industry (wood shred)
- Non-wovens (industrial fibres)
- Fibre composite industry (carbon fibres)

Anchor of wool application


Textile industry (wool)

- Applications: fineness measurement
- Parameters: fibre thickness



Wool standard sample.


Thickness Histogram



Histogram of wool thickness.
Anchor of cotton application


Chipboard and paper industry (wood shred)

- Applications: fineness measurement, length measurement
- Parameters: fibre thickness and length



Wood shred sample.


Counted particles : 35708
 min. diameter : 39.31 [µm]
 Median min. diameter : 24.26 [µm]
 Ellipticity : 2.3055, s: 1.0355
 Crystallinity : 1.1585, s: 0.0912
 Measured objects image fraction [%] : 18.52
 Object area fraction [%] : 62.42
 All objects image fraction [%] : 29.66



Percentile	Value [µm]
0.0	0.0
4.70	4.70
5.0	5.0
26.38	26.38
25.0	25.0
79.66	79.66
50.0	50.0
183.31	183.31
75.0	75.0
410.31	410.31
95.0	95.0
1074.98	1074.98
100.0	100.0
1521.77	1521.77

Hisam of wood shred thickness.

Counted particles : 35708
 max size : 76.80 [µm]
 Median max size : 48.60 [µm]
 area.wt max size : 294.53 [µm]
 Ellipticity : 2.3055, s: 1.0355
 Crystallinity : 1.1585, s: 0.0912
 Measured objects image fraction [%] : 18.52
 Object area fraction [%] : 62.42
 All objects image fraction [%] : 29.66




Percentile	Value [µm]
0.0	0.0
13.30	13.30
5.0	5.0
51.75	51.75
25.0	25.0
147.24	147.24
50.0	50.0
294.53	294.53
75.0	75.0
644.87	644.87
95.0	95.0
1604.35	1604.35
100.0	100.0
2481.74	2481.74

Histogram of wood shred length.
Anchor of fibre crimp application

Non-wovens (industrial fibres)


- Applications: fibre wave length measurement
- Parameters: crimp frequency (ASTM), crimp ratio (ASTM)



Fibre crimp sample.

parameters	mean value	standard deviation	total sums
avg. Wave length	0.37	0.10 cm	number of fibres:: 24
avg.fibre length:	3.62	0.40 cm	stretched fibre length 86.91 cm
projected length	2.79	0.24 cm	projected fibre length 66.93 cm
FERET length	2.79	0.24 cm	FERETMAX 67.08 cm
crimp degree	1.30		number of bows 377
crimp ratio	0.23 (ASTM)		
bow frequency	5.63	Bögen/cm	
crimp frequency	2.17 (ASTM)	1/cm	


Histogram of wave length distribution



Histogram of fibre crimp wave length.
Anchor of carbon fibres application


Fibre composite industry (carbon fibres)

- Applications: fineness measurement, length measurement
- Parameters: fibre thickness and length




Carbon fibres sample.

Counted objects : 328
 Fibre Thickness : 9.00, s: 2.40 [μm]
 Median Fibre Thickness : 8.73 [μm]
 lwt Fibre Thickness : 9.30 [μm]
 Measured objects image fraction [%] : 0.19
 All objects image fraction [%] : 3.68
 Curvature : 0.00056473
 Rectangularity : 0.9967, s: 0.0578
 Object area fraction [%] : 5.16



Histogram of carbon thickness.

Fibre length : 135.37 [μm]
 Median Fibre length : 118.67 [μm]
 Measured objects image fraction [%] : 0.19
 All objects image fraction [%] : 3.68
 Curvature : 0.00056473
 Rectangularity : 0.9967, s: 0.0578
 Object area fraction [%] : 5.16



Histogram of carbon length

Powdershape

1. The Powdershape system

Powdershape is a characterisation system for powders and generally all kinds of particles especially developed for quality inspection. It is unique in its set of features, simple user interface, and rapid characterisation process from the sample preparation to the printed certificate. The measurement range is from 1 μm to 20 μm particle size.

Two kind of Powdershape systems exist:

- Powdershape - where the samples are distributed stepwise by hand on the scanner Powdershape automatic - where the samples are feeded automatically by the sample feeder on the scanner. (sample feeder: patent pending).



The very simple handling of "Powdershape automatic"

1. Simply distribute your particle samples on the board
2. Enter a description and select grain size, max size or other reports with a single click!
3. All the specimen are scanned and characterized automatically. Obtain the final reports or insert results into quality chart database.

Your advantages:

- A very quick measurement method
- Easy and fast specimen/sample preparation
- Minimisation of possible operator errors
- Consistent and reproducible results
- Very attractive price/performance relation
- Reliable and rapid industrial method for quality inspection

he small brother of Powdershape "Microshape"

A new simply product "Microshape" has been issued for powder characterization that has just two measuring ranges

- o 2 - 150 μm using the digital microscope with 14600 dpi
- o 10 - 1000 μm using an 3600 dpi film scanner




2. Powdershape features

Key features

- Quick measurement of particle size and shape distribution
- Characterisation of optical properties: colour and transparency
- Ready-to-use report/certificate templates for grain size, mesh size and others.
- Quality charts
- Processing and storage of test data (for further external use and analysis)
- SAP compatibility available

Particles of a size range from 5 μm to 50000 μm (5 cm) can be measured.

Simple user interface



To facilitate the particle characterisation, only 3 steps are needed to obtain a certificate:

1. Image acquisition
2. Particle recognition and measurement
3. Analysis by means of statistical evaluation and report/certificate creation

Technical data

The Powdershape System consists of:

- The Powdershape software
- State-of-the-art image analysis computer running Microsoft Windows NT/2000/XP
- A slide and/or flatbed scanner certified for PowdershapeStandard office suite


3. Powdershape applications

- Sand
- WC-Co
- Nano particles
- High aspect ratio particles (gypsum)
- High aspect ratio particles (taken with a SEM)


Anchor of sand application

Sand

- Applications: particle sizer
- Parameters: grain size



Sand sample.




Sand grain size distribution.
Anchor of wc-co application

WC-Co

- Applications: particle sizer
- Parameters: grain size




WC-Co sample.



WC-Co grain size distribution.
Anchor of nano particles application


Nano particles

- Applications: particle sizer
- Parameters: grain size



Nano particles sample.

Counted particles : 425
grit size : 0.64, s: 0.30
Median grit size : 0.59
vol.wt grit size : 1.00
Ellipticity : 1.3106, s: 0.2254
Crystallinity : 1.0679, s: 0.0357




Percentile	[µm]
0.0	0.0
0.22	0.22
5.0	0.51
5.01	0.51
25.0	0.79
50.0	1.00
75.0	1.26
95.0	1.61
100.0	1.63

Size distribution histogram of nano particles.
Anchor of gypsum high aspect ratio particles application


High aspect ratio particles (gypsum)

- Applications: particle size, particle length
- Parameters: thickness, length, aspect ratio (calculated)



Gypsum sample.


Counted particles : 38115
min. diameter : 35.27 [µm]
Median min. diameter : 22.81 [µm]
Ellipticity : 2.2586, s: 1.0542
Crystallinity : 1.1477, s: 0.0843
Measured objects image fraction [%] : 19.29
All objects image fraction [%] : 25.10
Object area fraction [%] : 76.84



Percentile	[µm]
0.0	0.0
4.70	4.70
5.0	5.0
22.32	22.32
25.0	25.0
62.80	62.80
50.0	50.0
132.57	132.57
75.0	75.0
288.50	288.50
95.0	95.0
863.78	863.78
100.0	100.0
1443.18	1443.18

Histogram of gypsum thickness.

Counted particles : 38115
max size : 69.79 [µm]
Median max size : 43.37 [µm]
area.wt max size : 223.75 [µm]
Measured objects image fraction [%] : 19.29
All objects image fraction [%] : 25.10
Object area fraction [%] : 76.84




Percentile	[µm]
0.0	0.0
14.87	14.87
5.0	5.0
44.89	44.89
25.0	25.0
126.07	126.07
50.0	50.0
223.75	223.75
75.0	75.0
401.53	401.53
95.0	95.0
1154.41	1154.41
100.0	100.0
2081.39	2081.39

Histogram of gypsum length.
Anchor of sem high aspect ratio particles application


High aspect ratio particles taken with a SEM

- Applications: particle size, particle length
- Parameters: thickness, length, aspect ratio (calculated)




Stiks sample taken with a SEM.

Counted objects : 58
 Fibre Thickness : 1.86, s: 0.16 [µm]
 Median Fibre Thickness : 1.84 [µm]
 vol.wt Fibre Thickness : 1.88 [µm]
 Curvature : 0.00047027
 Rectangularity : 1.0061, s: 0.0396



Histogram of stick thickness.

Counted objects : 58
 Fibre length : 12.20, s: 1.89 [µm]
 Median Fibre length : 12.36 [µm]
 Curvature : 0.00047027
 Rectangularity : 1.0061, s: 0.0396



Histogram of stick length

Powdershape & Fibreshape

Ver. 2011

